

## BUFFER PREPARATION GUIDELINES FOR NANO-ITC

- Buffer with a low ionization enthalpy is best such as **PHOSPHATE, ACETATE or HEPES**.
- Buffers such as **TRIS** and **TRICENE** are not recommended if only one experiment is performed to ensure an accurate enthalpy.
- Usage of **SALTS** is not restricted, but higher salt concentrations may inhibit binding events if the ligand and macromolecule form ion pairs
- Any **pH above 4.0** is OK, however the pH of the cell and syringe must match within **0.1 pH** units
- **Reducing agents are OK**. **TCEP** is preferred to **DTT and B-mercaptoethanol** , but any of the three can be used if concentrations in the syringe and cell are matched **ABSOLUTELY**. The cell and syringe sample must be dialyzed at the same time in the exact same buffer, and the samples **MUST NOT** be degassed
- **DMSO** and **ETHANOL** are **OK** to use to solubilize ligands, however the DMSO or Ethanol concentration must be matched absolutely in the ITC cell as well, down to 0.01%. Extreme care must be taken when pipetting to ensure a matched system as the heat of dilution of these compounds is very large.
- **GLYCEROL, SUCROSE, UREA and DETERGENTS** at concentrations where they will not easily foam are **acceptable**.

Please remember that ITC is sensitive enough to measure the heat of dilution of most mismatched compounds at even very low concentrations. Concentrations in the cell and injecting syringe **MUST** be the same (IMIDAZOLE for example)

Compounds that are **NOT OK**, are **IODIDE** compounds and certain **CYANIDE** compounds as they leach gold, also acids at **pH lower than 4.0 CANNOT BE USED** as they are not compatible with the stainless steel syringe needle.